

No	Retrofit	Description	Category					Machine						Status	
			Performance	Quality	Security	Innovation	Cost saving	CL2	CL3	TLD	PAD	KAD	TPA		TC6i
1	Nonwoven Camera System	More reliable detection of the Nonwoven and Ribbon presence		x					x	x	x	x			available
4	Conversion Radialpress to Linearpress	Improving the stability of the tampon due to linear press type		x				x	x	x	x	x			available
5	Conversion tampon head colsing unit with folding link	To be able to reduce the wrapping foil width by using this new unit	x	x				x	x						available
14	CL3 Head forming unit 6 ways	Better cello head shape due to multiply head forming		x				x	x						available
15	Cello control camera system	More reliable cello detection with a camera system		x				x	x						available
16	String guide unit with ionizer	Better handling electrical-charged string		x				x	x	x	x	x			available
17	Perforation station	Perforation opening aid instead of teartape				x		x	x						available
19	Conversion complete head forming unit to CL3 design	More reliable unit for the head forming	x	x				x							available
24	Nonwoveneinheit with Vacuum	Better positioning of the Nonwoven on the ribbon		x				x	x	x	x	x			available
25	ATM	To prevent telescoping on the tampon		x	x	x		(x)	x	x	x	x			In progress
26	Upgrade localization Cello error	Fast detection on damaged insertion tubes		x				x	x						vorhanden